

In the operation of an autoclave when the discharge outlet valve is closed, ----- minutes are needed for sterilisation?

20

Which of the following equipment utilizes steam for sterilization?

Autoclave

The labels of stock solutions in the laboratory should contain all except one of these information

Name of person that prepared it

Dissecting trays contain-----

Wax

Used culture plates are better sterilized in the laboratory using-----

Oven

To display general instructions in the laboratory, which of the following factors should be least considered?

Wordings

In the laboratory disposal of spent materials are done on -----

The incinerator

Which of these items of the laboratory is out of place in the group?

Hot plate

The operational range of temperature in an incubator is between -----

22-50°C

The culture of microorganisms is best done in the laboratory using

Laminar flow

Whatsapp: 08089722160 or click here for TMA assistance

Practice E-exams & Chat with course mates on noungeeks.net