## Join group: T.me/NOUNSTUDENTSFORUM CLICK TO DOWNLOAD MORE TMA PQ

blowing any drop of the liquid to be transferred hanging on the top of the	e pipette.
is done to keep the interface between the lid and the body of desiccator air tight?	
greasing it	
Provided this apparatus, does not include those required in a set user ordinary filtration	
crucible tongs	
During titration it is advisable to hold the tap with your	_ hand
left	
The available sizes of volumetric flask is	
5cm3 to 5dm3	
The tap or stopcork of burette is lubricated with grease to prevent	
leakage	
Air tight container containing drying agent or desiccant is called	
desiccator	
The capacity of the flask are respectively indicated at thevolumetric flask	of the
neck of the flask	
The criteria for a precipitate to be dried ignited depend upon	
Nature of precipitate and type of filtering	
Volumetric flask are mostly used for preparing solution	
standard solution	
0/10	

Whatsapp: 08089722160 or click here for TMA assistance

Practice E-exams & Chat with course mates on noungeeks.net