Join group: T.me/NOUNSTUDENTSFORUM CLICK TO DOWNLOAD MORE TMA PQ

- 1. What is a lab coat?
- A coat made of white material worn in the lab
- 2. The essence of dissecting an animal or insect is
- To display the internal organs/system
- 3. Wearing of a lab coat is
- Professional
- 4. Wearing of a lab coat is for
- All of the above
- 5. Animals for dissection are held down on a board or tray with

Pins

6. _____constitute the largest number of living vertebrates

The bony fish

7. To maximize the light gathering ability of the microscope lens, you add_____between the object and the coverglass

Immersion oil

8. Why is wearing of short is the lab not advisable?

To avoid injury

9. You open up a rat with

Scissors

10. Which of the following statements is incorrect about a microscope?

It is kept in a damp place

Whatsapp: 08089722160 or click here for TMA assistance

Practice E-exams & Chat with course mates on noungeeks.net